

## MODALIZER+ V6 Release Notes

Version 6.0.1.0, Sep 10, 2021

# **Summary**

A complete new graphic design with major improvements and fixes. Full support for Secure DICOM networking. DICOM templates for every media type. Subscription model.

V6.0.1.0 Sep 10, 2021 Build 479

MPLUS-197 Shorten application start time by deferring calls to license web service when the license is still valid

V6.0.0 Aug 31, 2021 Build 473

## Bug

MPLUS-196 Your license will expire in 29 days message shows for subscription license.

V6.0.0 Jul 27, 2021 Build 472

MPLUS-193 DICOM Viewer not working - missing resource file

MPLUS-192 When running DICOM Viewer only don't check activation WS

MPLUS-190 if activation web service cause error the modalizer crash

V6.0.0 Jul 22, 2021 Build 470

New features and improvements



MPLUS-172 Uninstall options data options. Uninstall now lets the user decide weather to keep or discard the data.

MPLUS-138 Update to .net version 4.5. Supported OS - Windows 10 and newer.

MPLUS-132 Update security settings.

MPLUS-130 Re-order main screens button

MPLUS-129 Improvements on Settings screen

MPLUS-117 Add user manual to the setup

MPLUS-114 Local Archive screen improvements

MPLUS-112 License update flow. License is now in subscription model with automatic online renewal.

MPLUS-109 Add DICOM C-GET.

MPLUS-108 Add option to select certificate in CurrentUser OR LocalMachine store by user selection

MPLUS-107 Add Import files action to Local Archive Screen

MPLUS-101 Open DCXSEC.MutualAuthentication properties which is now enabled with last RZDCX build

MPLUS-100 Add Listener indication on main screen

MPLUS-96 Remote DICOM AE Edit popup security options

MPLUS-95 DICOM Listener security settings

MPLUS-94 Add certificate browsing option in settings

MPLUS-88 Multiple templates for dicom files in dicomizer based object class

MPLUS-5 Add RZDCX version to home screen info tooltip

MPLUS-181 Remove taken images and files on DONE (and not on App close)



MPLUS-147 Add DCXAPP.ProjectVersion to tooltip in home screen so it shows exact RZDCX build

MPLUS-125 MODALIZER 6 New UI/UX

MPLUS-106 Create x64 Version

MPLUS-184 Evaluation form Name field is Mandatory

MPLUS-174 Add progress bar to export

### Bug fixes

MPLUS-180 Last taken file is not deleted

MPLUS-179 Q/R Doesn't fill local archive if local archive was not enabled at application start

MPLUS-177 Files from Q/R are stored in LocalStorage folder and never deleted

MPLUS-173 Desktop shortcut icon doesn't work

MPLUS-171 Secure MPPS

MPLUS-166 query retrieve - search doesn't work with patient name, works only with patient id

MPLUS-165 Click cancel after clicking Browse on any certificate select shows error

MPLUS-163 Fix security configuration - loading and saving

MPLUS-162 MODALIZER+ leaves empty study folders in ProgramData

MPLUS-161 MODALIZER+ db file is left in program data after uninstall

MPLUS-160 CD Viewer doesn't work

MPLUS-159 C-GET fails when use UNICODE is on

MPLUS-158 C-GET bug: can't log wide string filename

MPLUS-157 Bottom strip invisible in the contrast presets list control

MPLUS-155 Conversion fails for this PDF

MPLUS-154 Radio buttons state on Convert screen is not saved between sessions



MPLUS-151 MR Image Sorting within frame fails on image position patient sign change

MPLUS-149 M+ v6 new UI issues (duplicate)

MPLUS-146 PDF extract in x64 bit fails

MPLUS-144 Tab order in new/edit remote server not stopping in Retrieve method

MPLUS-143 Certificate selection UI always creates new certificates instead of using the one selected by the user

MPLUS-142 X64 version installed in Program Files (x86)

MPLUS-141 Verification method get into effect only after save (5.5 and 6.0)

MPLUS-140 MODALIZER 5.5 - Echo result indicator is behind cancel button in Edit Storage Server

MPLUS-137 Secure listener doesn't work

MPLUS-136 Add button of thumbprints should give text box and not open browse store

MPLUS-135 Listener echo button fails the security is on

MPLUS-131 Typo in tooltip of New Patient button

MPLUS-127 Crash on Mouse down (left button) in Local Archive screen

MPLUS-126 First run after install gives this error. Next run works.

MPLUS-124 Crash in Local Archive screen

MPLUS-122 Two step verification in MWL

MPLUS-121 "Config per user" is not stored in the MODALIZER+

MPLUS-120 MODALIZER release build doesn't have digital signature

MPLUS-119 Crash on open Local Archive if DB is empty

MPLUS-113 Unrecognized attribute on downgrade

MPLUS-111 RZDCX version installed with MODALIZER ad MODALIZER+ is wrong



MPLUS-103 On first opening of Settings screen "Expanded" icon is not shown on Storage servers and Worklist servers captions

MPLUS-99 TLS control are enabled on Setting form opening even if "Use security" is not checked

MPLUS-97 Listener gets stuck with 3 repeating echos

MPLUS-90 Show header does not display empty sequences

MPLUS-89 Incorrect error message when Q/R service is not supported by PACS

MPLUS-81 Edit DICOM file and save to new folder doesnt work

MPLUS-44 Change Buy Now link

MPLUS-42 Message boxes jump behind main window